


Product Data Sheet

PIA-PINK-Block



Blocking solution optimized to prepare blot membranes for PIA-PINK immunostaining.

Products

Type	Name	Icon	Size	Product Number
ready-to-use reagent	PIA-PINK-Block		5 vials 10 mL each	# 100-01

Overview


Name	PIA-PINK-Block
Description	Ready-to-use reagent for PIA-PINK immunostaining. PIA-PINK-Block contains serum albumin, buffer, salt and additives in a composition optimized to prepare blot membranes for immunostaining with PIA-PINK solution.
Tested Applications	Western Blot, Dot-Blot
Other Applications	ELISA, electron microscopy
Product Size	Box containing 5 vials 10 mL ready-to-use reagent, each. 10 mL reagent is the recommended volume to cover a Midi or Mini Blot membrane.
Delivery	At ambient temperature (8 °C - 30 °C)
Storage	At room temperature (15 °C - 25 °C) or refrigerated (4 °C - 8 °C). Do not freeze. Ideally the vials stay horizontally and protected from sun light.
Shelf life	6 months. Use-up immediately when a vial is opened.

For in vitro research use only. Not for diagnostic or therapeutic use.

PIA-PINK are reagents used in immunological applications without any further dilution.

PIA-PINK-Block is used for immunostaining with PIA-PINK reagents, like PIA-PINK-Mouse, -Rabbit or -His. Please refer to the respective protocols for further instructions.

Valid from 2021-10-25



PiNa-Tec®

Notkestraße 85 - D-22607 Hamburg, Germany

Telephone: +49 (0) 40 646 33960 - Mobil: +49 (0) 176 209 40402

Email: katja.werner@pina-tec.de - Internet: www.pina-tec.de